

Intro to spatial omics

04/22/2025

Jiwoon Park, Ph.D.

We want the whole platter and their interactions

Bulk



Counts/group

Single cell



Counts/cell



Spatial



Counts/cell + location



Before spatial “omics”

Major events in evolution of prequel techniques




Molecular biology with in-situ hybridization

Major prequel WM ISH atlases




Sequencing and multiplexing

Major events in evolution of current-era techniques



Sequencing and multiplexing



Explosion since 2020

The image displays two screenshots of the Nature Methods journal website. The top screenshot shows the article "Method of the Year: spatially resolved transcriptomics" by Vivien Marx, published on 06 January 2021. It has 170k accesses and 454 citations. A "Publisher Correction" notice is visible. The bottom screenshot shows the article "Method of the Year 2024: spatial proteomics" by Vivien Marx, published on 06 December 2024. It has 62k accesses and 420 altmetric mentions. A "Publisher Correction" notice is also visible. Both articles are part of the "Method of the Year" series.

nature methods

Explore content ▾

About the journal ▾

Publish with us ▾

nature > nature methods > editorials > article

Editorial | Published: 06 December 2024

Method of the Year 2024: spatial proteomics

Save

Related Papers

Chat with paper

Nature Methods 21, 2195–2196 (2024) | Cite this article


62k Accesses | 420 Altmetric | Metrics

Approaches for profiling the spatial proteome in tissues are the basis of atlas-scale projects that are delivering on their promise for understanding biological complexity in health and disease.




Humans by nature are inquisitive. We love to explore, and as we do, we make maps. Perhaps it is then not surprising that we find ourselves in an age of atlas-building as we explore the incredible complexity of biological systems with cutting-edge methods to better understand how cells, tissues, organs and organisms connect structure and function.




Landscape of Spatial Technologies




Every organ has a spatial organization



Spatial organization of cell types and expressions matter




Liver fibrosis development




Wound healing and scar formation

Spatial omics created opportunities for digital pathology

19th Century Histopathology




Tissue morphology,
TNM staging/grading




One treatment for all


21st Century Molecular Pathology



Companion diagnosis



NOW & FUTURE Digital Pathology





Spatial + Multi omics




Personalized diagnostics

Crash course to spatial omics for today




If you want to run spatial omics...


Question: Why do you need spatial omics for your project?




Before starting your spatial journey





Spatial omics workflow, simplified



Factors to consider every step



Considerations during sample preparation: fixation



Different protocols but check:

- Uniform penetration,
- Temperature,
- Fixative concentration & time


pH matters!

Considerations during sample preparation: preservation

Table 1	 FFPE Tissue	 Frozen sections
<i>Tissue morphology preservation</i>	✓ Yes	✗ Poor
<i>Long term storage</i>	✓ Several years at RT	✗ 1 year at -80°C
<i>Section thickness</i>	✓ Thin, 1-50 µm	✗ Thicker → lower resolution
<i>Enzymes activity and labile antigens preservation</i>	✗ No	✓ Yes
<i>Antigen immunoreactivity</i>	✗ Epitope masking, antigen retrieval often necessary	✓ No epitope masking
<i>Fixation</i>	✓ Prior to embedding	✓ After sectioning

- Antigen retrieval protocol may vary
- Decalcification protocol
- Tissue clearing for some platforms
- Some tissues tend to be more sensitive to degradation (enzymatic activities)

Factors to consider every step



Available platforms to consider



CosMx (EIPM)



IMC (EIPM)



CellScope



CODEX
(Phenocycler)



Xenium



VisiumSD/HD



Stereo-seq




PIXELGEN
TECHNOLOGIES

Questions?

Imaging-based spatial transcriptomics


Imaging based spatial transcriptomics



• ***In situ* sequencing (ISS)** involves hybridizing probes to target transcripts, amplifying the signal through techniques like rolling circle amplification, and then sequentially reading nucleotide sequences using iterative rounds of fluorescent labeling and imaging.

• ***In situ* hybridization (ISH)**, detects specific RNA or DNA sequences by hybridizing many fluorescently labeled probes directly to target molecules, enabling spatial localization.

Basic concept for fluorescence in situ hybridization (FISH)



Challenges:


1. **Multiplexing constraints:** Traditional fluorescence microscopes can only detect a few fluorophores at once, limiting researchers to visualizing only a **handful of RNA species** in a single experiment. What if scientists wanted to track **thousands of genes simultaneously**?
2. **Optical crowding:** In densely packed cellular environments, RNA molecules can cluster so tightly that **fluorescence signals overlap**, making it difficult to distinguish and count individual molecules.

Barcoding and detection approach vary by technology



Signal amplification

Error correction

Plexity increase




In Situ Sequencing (ISS)



low gene detection efficiency
(low efficiency of *in situ* cDNA synthesis)


Improvements in the workflow

Expansion sequencing (ExSeq)




Hydrogel to solve optical crowding (Chemical anchor and swelling the tissue) -> increased detection efficiency.

STARmap – SNAIL probe




No RT, better capture efficiency
Later STARmap PLUS & RIBOmap with various probe design optimization efforts.

FISSEQ




Sequencing by ligation;
untargeted, high multiplexing capacity but low efficiency.

Going 3D – coming soon




STARmap: discovery and distribution of cell types in 3D



Imaging-based spatial proteomics

Immunostaining to detect protein in the lab




- Spectral overlap between fluorophores, limiting the number of detectable targets
- Cross-reactivity of antibodies, leading to non-specific staining
- Limited multiplexing capacity, restricting simultaneous detection of multiple markers
- Tissue autofluorescence, reducing signal clarity in fluorescence-based IHC


Immediate solution: remove the tags after use

Cyclic Multiplexed Immunofluorescence Methods

Fluorescent cyclic approaches




Fluorophore labelled (indirect immunofluorescence)



Signal removal with:


- fluorophore inactivation by oxidation (e.g., hydrogen peroxide)
- Alkaline pH
- Chemical bleaching (e.g., LiBH_4 /lithium borohydride)
- Light (photobleaching)
- Antibody elution buffers

Adding DNA reporter tags to improve detection



- After this initial staining, the detection process uses multiple rounds of imaging without additional antibody applications.
- Each imaging cycle, fluorescently labeled reporter probes complementary to specific DNA tags hybridize to their target sequences, generating detectable signals.
- When these hybridization events occur, the resulting fluorescent patterns serve as unique identifiers for different proteins.
- Between imaging cycles, the reporter probes are removed through rapid chemical or enzymatic processes, allowing new probes to bind in subsequent rounds.


Mass Cytometry readout instead of imaging




- **High Multiplexing:** protein targets (up to 40 or more) in a single tissue section.
- **Preservation of Tissue Architecture:** IMC preserves the spatial organization and relationships of cells within the tissue, unlike flow cytometry which requires cell dissociation.

- **Lower Resolution:** relatively lower (1 μ m) compared to other proteomics imaging techniques
- **Limited Availability of Isotopes:** The number of available isotopes for staining can be a limiting factor in expanding the number of markers that can be simultaneously analyzed. In some cases, metal contamination from therapeutic or diagnostic agents can interfere with IMC signal detection.
- **Cost:** Reagents (e.g., metal-tagged antibodies) can be expensive, and difficult to optimize (and harmful for the environment)


Same concept for metabolomics



Spatial proteomics development for better detection



Spatial proteomics



Going for multi-omics

Spatial Proteomics: Multiomics Capabilities



Manual Sequential Workflows



Phenocycler-Fusion
Akoya Biosciences



Phenomager
Akoya Biosciences



Cell DIVE
Leica Microsystems



CellScape
Canopy

Automated RNA & Protein Co-detection



COMET
Lunaphore



MACSima
Miltenyi Biotec




Hyperion Imaging Systems
Standard BioTools Inc.




InSituPlex
Ultivue (merged with Vizgen)

Protein and Other Metabolites Profiling




MALDI HPLC-IMS
Bruker

Spatial metabolomics + proteomics



And beyond traditional proteomics: Molecular pixelation


a



- 1. Protein abundance:** Similar to flow cytometry or CyTOF, measuring the overall expression levels of proteins per cell.
- 2. Individual protein spatial distribution:** Quantifying whether a single protein is randomly distributed, polarized, or clustered within the cell.
 - Polarity scores indicate whether a protein is organized (<0), clustered (>0), or randomly dispersed ($=0$).
 - Example: T cells stimulated via CD3 show polarized CD3 localization on one end of the cell.
- 3. Protein colocalization analysis:** Determining how different proteins are positioned relative to each other on the same cell.
 1. Colocalization scores range from <0 (segregated proteins) to >0 (proteins physically co-localized).
 2. Example: Rituximab treatment of RAJI cancer cells polarizes CD20, a mechanism critical for CAR-T therapy and antibody-based drugs.

Sequencing-based spatial transcriptomics


Sequencing based spatial transcriptomics




- Method that ultimately uses sequencing to decode spatial information: Involves library preparation & sequencing
- Spatial information is in the barcodes, molecular information is in the sequencing library
- Varied by **capture methods (spatial indexing and/or ligation of spatial barcodes)**
 - Capture on barcoded surface
 - Applying barcode via device
 - Apply barcode directly on tissue

Capture method varies within sequencing-based spatial omics

Spatial transcriptomics (sequencing)



First version of spatial sequencing: 2016, Ståhl et al.




- ~200 million oligonucleotides across 1007 features
- 100 mm diameter, 200 mm center-to-center distance
- capture area of 6.2 X 6.6 mm² (100-200 µm resolution)

Capture rate were much lower than mFISH methods

Each Oligonucleotide on barcoded slide contained:

1. A **cleavage site** to release captured molecules during processing,
2. An **amplification and sequencing handle** for downstream library preparation,
3. A **spatial barcode** encoding positional information,
4. A **unique molecular identifier (UMI)** to track individual transcripts and minimize PCR biases, and
5. An **mRNA capture region** with a poly(dT) tail to hybridize polyadenylated mRNA from the tissue.


Visium (2019), first commercialization of the technology



- resolution that is **2-4X higher**: The slide contains approximately **5,000 spots**, each with a **55 μm** diameter
- Offering FFPE support, streamline with CytAssist (2022)



Visium workflow with CytAssist



Facilitate transfer of transcriptomic probes in FFPE or fresh frozen samples with Visium CytAssist. In the Visium CytAssist workflow, sectioning, tissue preparation, staining (H&E or IF), and imaging take place on a standard glass slide. After probe hybridization, two standard glass slides and a Visium slide with two Capture Areas are placed in the CytAssist instrument so that the tissue sections on the standard slides can be aligned on top of the Capture Areas. Within the instrument, a brightfield image is captured to provide spatial orientation for data analysis, followed by hybridization of transcriptomic probes to the Visium slide. The remaining steps, starting with probe extension, follow the standard Visium workflow outside of the instrument.

Slide-seq commercialization with Curio for higher resolution

Slide-seq (Rodrigues et al.) increased capture efficiency by applying 10 μm beads onto a thin glass

Tissue section is placed onto
spatially indexed beads and
RNA is captured

cDNA
synthesis

Bead recovery
from substrate

NGS library
preparation

Sequencing


Map of differential
expression

3 mm x 3 mm


3 mm x 3 mm

2024: True single-cell resolution after optimization (Curio Trekker) by directly diffusing barcodes to cell nuclei using UV


- + No instrument needed, simple workflow
- Sequencing cost, only FF





Pushing for high resolution with high capture rate



Similar approach with VisiumHD




two $6.5 \times 6.5 \text{ mm}$ capture areas within an $8 \times 8 \text{ mm}$ fiducial frame, densely packed with **~11 million uniquely barcoded $2 \times 2 \mu\text{m}$ squares** (typically analyzed at $8 \times 8 \mu\text{m}$ bins)




More sub-micron spatial transcriptomics on the horizon

Stereo-seq (STOmics, 220 nm)



Illumina Spatial (Seq-Scope, 0.5 - 2 μm)





Now with better sequencing chemistries

Element AVITI24



Avidites leverage the power of multivalent binding to reduce costs while increasing accuracy




Singular Genomics G4X





Multomics
RNA, Protein, fH&E

High-Throughput
20x More Samples per Week

Fast
3-Days Sample-to-Discovery




Multi-omics sequencing – DBiT-seq




Spatial barcoding with photo-activation


GeoMx DSP




Light-Seq




Temporal measurements



Perturb-Map




Perturb-map




Questions?

Factors to consider every step




Spatial omics data analysis

Data processing + analysis



combine as cell metadata for mapping
(single cell masks/expression profiles)

Cell type and expression profiling + Tissue microenvironment characterizations




Typical single cell
transcriptome analysis




cellular neighborhood/interaction analysis

Cell segmentation – why do we care?





- **Understanding Cellular Interactions:** Cell segmentation allows scientists to explore cellular interactions in normal and pathological conditions, enhancing our understanding of fundamental biological processes.
- **Analyzing Biological Features:** Cell segmentation empowers scientists to analyze essential biological features like cell count, type, shape, and more. This analysis provides insights into how these features change over time and in response to different conditions.
- **Capturing Biologically Relevant Information:** The morphology of cells reflects their physiological state, and effective cell segmentation captures biologically relevant morphological information.
- **Accelerating Drug Discovery:** Insights gained from cell segmentation can drive drug discovery, diagnostics, and various other critical areas in biology, pharmacology, and personalized medicine.

Challenges in cell segmentation




- **Heterogeneous Shapes:** Cells come in diverse and dynamic shapes, making it nearly impossible to define mathematical shape models.
- **Variation in Size and Shape:** Unlike nuclei, which are typically similar-sized structures, the cytoplasm exhibits significant variations in shape and size.
- **Weak Boundary Gradients:** Cells that are in close proximity can have weak boundary gradients, making it challenging to separate them.
- **Makeshift Nature of Segmentation Approaches:** The applicability of segmentation methods can be limited by dataset constraints, including differences in staining or imaging modality, artifacts in image capture, or morphological differences.

Segmentation-free methods if interaction is not important



Considerations for spatial data analysis workflow



QC parameters (CosMx)

Protein

Gating and comp. with neg probes

1. Neg probes

- Keep cells with neg probes between 2 and 50 (too low is indication of data preprocessing issue, too high means high background -> unreliable intensity values)

2. High expression

- Remove cells where at least X% of proteins are in the Yth percentile (default X = 50, Y = 90). Likely autofluorescence cells

3. Low expression

- Remove cells where fewer than V proteins are in at least Wth percentile (default V = 3, W = 50). Likely not real cells.

RNA

Data modeling and filtering

1. Cell QC

- Min counts per cell (20)
- Ratio of negative to targets counts (0.1)
- Counts per unique gene (“complexity”) (1)
- Removal of unusually large or small cells (p value of Grubb’s test = 0.01)

2. Field of View (FOV) QC

Removing low-quality FOVs: (based on mean)

- low average counts per cell (10)
- counts too close to background (0)


3. Target QC

- Score test thresholds for detecting genes above background (0.01 is the score test p-value for detection over background)
- Plot and confirm (I’m checking what AtoMx’s standard QC plot is)

Spatial analysis and data modeling


A

Spatial pattern discovery




B

RNA kinetics modeling




C

Multimodal integration



Need for large-scale datasets in production for reference



**HUMAN
CELL
ATLAS**

